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Exercise Accompanying the MS Tutorial

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Fig. 1: Persulfonated first-generation dendrimer with TREN core

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Fig. 3: Structure-perfect (top left) persulfonylated G1-POPAM dendrimer and defects.

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Fig. 4: MS/MS mass spectra of the first four fractions of the dendrimer synthesis. For the ease of analysis, some of the mass differences are given.

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Fig. 1 shows a small dendrimer persulfonylated at the periphery which is based on a TREN core ($C_{48}H_{54}N_4O_{18}S_6$) and an analogue ($C_{39}H_{43}N_3O_{15}S_5$), which lacks one branch.

- a) With which ionization method would you try to ionize the molecules? Provide reasons for your decision! Which ions do you expect to observe?

Fig. 2 shows an ESI-FTICR mass spectrum (top; spray solvent: MeOH) and below a MALDI-TOF mass spectrum (matrix: DHB), both from the exact same sample. As you may easily see, the two spectra differ much. In the MALDI mass spectrum, significantly more signals appear as compared to the ESI mass spectrum.

- b) If you take a closer look at the isotope pattern of the base peak, you may see that the third peak in agreement with the calculated pattern is quite intense. Explain why!
- c) Analyze the differences between ESI and MALDI mass spectra and draw conclusions with respect to the purity of the dendrimer! Which impurities are in the sample?
- d) Try to explain the additional signals observed in the MALDI mass spectrum! How could you distinguish fragments originating from gas-phase fragmentations from decomposition reactions occurring during the MALDI process and from impurities in your sample?

In Fig. 3, you see a series of POPAM dendrimers (1st generation) with different defects (marked gray). The sulfonyl groups were attached to the POPAM scaffold in a final reaction step by reacting the amino-terminated POPAM dendrimer with tosyl chloride under basic conditions. During chromatographic purification, you find several product fractions. Since all dendrimers contain repetitively the same structural elements, an NMR analysis remains ambiguous. You do not obtain secure information in which fraction which defect is present. You decide to use ESI-FTICR mass spectrometry and to obtain structural information by colliding protonated dendrimer ions in a collision experiment and analyzing their fragmentation behavior. Fig. 4 shows the result for the the first four fractions.

- e) Which MS/MS mass spectrum corresponds to the structure-perfect dendrimer? Suggest a plausible fragmentation mechanism which rationalizes the observed fragmentation reaction!
- f) If you consider the same fragmentation mechanism to occur for the other defects, how can you assign the other spectra to different defects? In the bottom spectrum, one observes many other signals as well. The complexity is a hint: At what? Consider which defect could be present in the last of the four fractions and provide reasons for your assignment!

Literatur/Literature:

- 1 *A Combined ESI- and MALDI-MS(/MS) Study of Dendrimers Persulfonated at their Periphery: Falsely Negative Results by MALDI MS and Analysis of Defects*, T. Felder, C.A. Schalley, H. Fakhrnabavi, O. Lukin, *Chem. Eur. J.* **2005**, *11*, 5625.
- 2 *How Useful is Mass Spectrometry for the Characterization of Dendrimers? "Fake Defects" in the ESI and MALDI Mass Spectra of Dendritic Compounds*, B. Baytekin, N. Werner, F. Luppertz, M. Engeser, J. Brüggemann, S. Bitter, R. Henkel, T. Felder, C.A. Schalley, *Int. J. Mass Spectrom.* **2006**, *249-250*, 138.
- 3 *Designer Dendrimers: Oligosulfonimide Dendrimers and Dendrons with Controllable Molecular Architecture*, O. Lukin, V. Gramlich, R. Kandre, I. Zhun, T. Felder, C.A. Schalley, G. Dolgonos, *J. Am. Chem. Soc.* **2006**, *128*, 8964.
- 4 *Mass Spectrometry as a Tool in Dendrimer Chemistry: From Self-Assembling Dendrimers to Dendrimer Gas-Phase Host-Guest Chemistry*, C.A. Schalley, B. Baytekin, H.T. Baytekin, M. Engeser, T. Felder, A. Rang, *J. Phys. Org. Chem.* **2006**, *19*, 479.